


















881VR Series Continuous / Intermittent				882VR Series Continuous-only			
Clockwise							
	881VR-160-1L pediatric	881VR-200-1L adult	881VR-300-1L adult	882VR-160-1L pediatric	882VR-200-1L adult	882VR-300-1L adult	882VR-760-1L surgical
							
Counter-Clockwise	881VR-160-1.5L pediatric	881VR-300-1.5L adult	882VR-160-1.5L adult	882VR-300-1.5L adult	882VR-760-1.5L surgical		
Digital Gauge							
	881VR-160-2.5L pediatric	881VR-300-2.5L adult	882VR-160-2.5L pediatric	882VR-300-2.5L adult	882VR-760-2.5L surgical		
	Mode Selection:						
OFF	√	√	√	√	√	√	√
CONT	√	√	√	√	√	√	√
INT	√	√	√	X	X	X	X
FULL	X	X	X	X	√	√	√
Regulated Vacuum Range	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~760 mmHg
Gauge Display Range	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~760 mmHg
Gauge Accuracy	± 4%	± 4%	± 4%	± 4%	± 4%	± 4%	± 4%
Intermittent Time	Suction: 16 seconds, Stop: 8 seconds			/	/	/	/
Flow Rate	>50 l/min (CONT) 7-9 l/min (INT)	>55 l/min (CONT) 7-9 l/min (INT)	>60 l/min (CONT) 7-9 l/min (INT)	>50 l/min	>55 l/min	>60 l/min	>70 l/min
Weight	1 lb / 0.455 kg			0.81 lb / 0.37 kg			
Dimensions	6.3" H × 3.15" W × 3.94" D / 16 × 8 × 10 cm			6.3" H × 3.15" W × 3.94" D / 16 × 8 × 10 cm			